

SGS Canada Inc.

P.O. Box 4300 - 185 Concession St. Lakefield - Ontario - KOL 2HO

Phone: 705-652-2000 FAX: 705-652-6365

Works #: 500286079 Project #: SFIS #:14365

03-May-2016

PWBS-Brant Haldimand Norfolk DCSB (Reg243) St. Basils

Attn: Paul Viszy (PWBS)

535 Paris Rd.. Paris, ON N3L 3E1,

Phone: 519-753-3266 Fax:519-753-7072

Date Rec.: 27 April 2016

LR Report: CA15591-APR16

Reference: co-lcated w/ Walter Gretzky

Copy: #1

CERTIFICATE OF ANALYSIS **Final Report**

Sample ID	Sample Date & Time	Temperature upon Delivery @ London Lab °C	Lead ug/L
1: Analysis Start Date			29-Apr-16
2: Analysis Start Time			18:37
3: Analysis Approval Date			02-May-16
4: Analysis Approval Time			10:12
5: MAC			10
6: MDL			0.01
7: TAP Standing - St.B-side Room 2030	23-Apr-16 08:15	14.3	0.23
8: TAP Flushed - St.B-side Room 2030	23-Apr-16 08:50	14.3	0.19
9: TAP Standing - WG-side Room 1018	23-Apr-16 08:30	14.3	0.39
10: TAP Flushed - WG-side Room 1018	23-Apr-16 09:05	14.3	0.15
11: TAP Standing - WG-side 2055 Fountain	23-Apr-16 08:14	14.3	0.02
12: TAP Flushed - WG-side 2055 Fountain	23-Apr-16 08:49	14.3	0.03
13: TAP Standing - St.B-side 1037 Fountain	23-Apr-16 08:30	14.3	0.05
14: TAP Flushed - St.B-side 1037 Fountain	23-Apr-16 09:05	14.3	0.02
15: TAP Standing - St.B-side 2056 Fountain	23-Apr-16 08:21	14.3	0.01 <mdl< td=""></mdl<>
16: TAP Flushed - St.B-side 2056 Fountain	23-Apr-16 08:56	14.3	0.02
17: TAP Standing - WG-side 1021 Fountain	23-Apr-16 08:32	14.3	0.01 <mdl< td=""></mdl<>
18: TAP Flushed - WG-side 1021 Fountain	23-Apr-16 09:07	14.3	0.02

MAC - Maximum Acceptable Concentration MDL - SGS Method Detection Limit

Patti Stark

Project Specialist Environmental Services,

Analytical



SGS Canada Inc.

P.O. Box 4300 - 185 Concession St. Lakefield - Ontario - KOL 2HO

Phone: 705-652-2000 FAX: 705-652-6365

Works #: 500286079 LR Report : CA15591-APR10 Project #: SFIS #:14365

Method Descriptions

Parameter	Description	SGS Method Code	Reference Method Code
Lead	Lead by ICP-MS Drinking Water	ME-CA-[ENV]SPE-LAK-AN-006	SM 3030/EPA 200.8